



SAGES

Laparoscopy Troubleshooting Guide

Society of American Gastrointestinal Endoscopic Surgeons

Developed and Distributed by the SAGES Continuing Education Committee

To minimize equipment malfunction, scheduled routine maintenance should be in place for all components of laparoscopy. Manufacturers' recommendations for routine replacement of some parts (e.g. bulbs) should be taken into consideration.

PREOPERATIVE PRECAUTIONS

Circulator Duties or Tasks

1. Assure table tilt mechanism is functional. Table & joints level, kidney rest down.
2. Consider using foot board and extra safety strap.
3. Position patient properly on O.R. table for cholangiography.
4. Assure notification of radiology technologist with time estimate.
5. Assure proper mixing and dilution of cholangiogram contrast solution for adequate image. For surgeons utilizing fluoroscopy for cholangiography, the patient should be placed on a table capable of supporting this task and appropriate shielding should be available.
6. Assure availability of Foley catheter and N/G tube.
7. Assure all power sources are connected and appropriate units are switched "on" (Don't use multi-socket single source or the circuit will overload).
8. Assure adequate volume of compressed gas (at insufflator and pressure irrigator). Backup full tank must be available. Ensure wrench and gasket are available.
9. Assure insufflator alarm is set appropriately. Assure tight connection between insufflator tubing and Luer-lock adapter.
10. Assure full volume in irrigation fluid container ...

11. (recheck during case). Check the electrosurgical unit; make sure auditory alarm of machine is functioning properly and the grounding pad is appropriate for the patient.
12. Check Veress needle for proper plunger/spring action and assure easy flushing through stopcock and/or needle channel.
13. Assure closed stopcocks on all ports.
14. If utilizing the gasless technique, assure that the operating room table has side arms capable of supporting the abdominal lift unit and that appropriate blades for the unit are available.

Scrub Tech/ RN Duties

15. Check sealing caps for cracked rubber, stretched openings.
16. Check to assure instrument cleaning channel screwcaps are in place.
17. Assure free movement of instrument handles and jaws.
18. If Hasson cannula to be used, assure availability of stay sutures and retractors. Check valves, plunger, spring, and seals on reusable Hasson cannulae.
19. Assure adequate printer film and video tape if documentation is desired.
20. Periodically send scissors and reusable trocars for sharpening.

PROBLEM

CAUSE

SOLUTION

1. Poor Insufflation/loss of pneumoperitoneum	CO2 tank empty Accessory port stopcock(s) not properly adjusted Leak in sealing cap or stopcock Excessive suctioning Loose connection of insufflator tubing at source or at port Hasson stay sutures loose Tubing disconnection from insufflator Flow rate set too low	Change tank Inspect all accessory ports. Open or close stopcock(s) as needed Change cap or cannula Allow time to reinsufflate Tighten connections Replace or secure sutures Connect tubing Adjust flow rate
2. Excessive pressure required for insufflation (initial or subsequent)	Veress needle or cannula tip not in free peritoneal cavity Occlusion of tubing (kinking, table joints, etc.) Port stopcock turned off Patient is "light" Cannula tip not in peritoneal space	Reinsert needle or cannula Inspect full length of tubing. Replace with proper size as necessary Fully open stopcock Give more muscle relaxant Advance cannula under visual control
3. Inadequate lighting (partial/complete loss)	Loose connection at source or scope Light is on "manual-minimum" Bulb is burned out Fiber optics are damaged Automatic iris adjusting to bright reflection from instrument Monitor brightness turned down Room brightness floods monitors	Adjust connector Go to "automatic" Replace bulb Replace light cable Re-position instruments, or switch to "manual" Readjust setting Dim room lights

SAGES thanks Valleylab, Inc. for their generous educational grant in support of the development and printing of this chart.

PROBLEM	CAUSE	SOLUTION
<p>4. Lighting too bright</p>	<p>Light is on "manual-maximum"</p> <p>"Boost" on light source is activated</p> <p>Monitor brightness turned up</p>	<p>Go to "automatic"</p> <p>Deactivate "boost"</p> <p>Readjust setting</p>
<p>5. No picture on monitor(s)</p>	<p>Camera control or other components (V.C.R., printer, light source, monitor) not "on"</p> <p>Cable connector between camera control unit and/or monitors not attached properly</p> <p>Cable between monitors not connected</p> <p>Input select button on monitor doesn't match "video in" choice</p>	<p>Make sure all power sources are plugged in and turned on</p> <p>Cable should run from "video out" on camera control unit to "video in" on primary monitor. Use compatible cables for camera unit and light source.</p> <p>Cable should run from "video out" on primary monitor to "video in" on secondary monitor</p> <p>Assure matching selections</p>
<p>6. Poor quality picture</p> <p>a. fogging/haze</p> <p>b. flickering, electrical interference</p> <p>c. blurring, distortion</p>	<p>Condensation on lens from cold scope entering warm abdomen</p> <p>Condensation on scope eyepiece, camera lens, coupler lens</p> <p>Moisture in camera cable connecting plug</p> <p>Poor cable shielding</p> <p>Insecure connection of video cable between monitors</p> <p>Incorrect focus</p> <p>Cracked lens, internal moisture</p> <p>Too grainy</p>	<p>Gently wipe lens on viscera; use anti-fog solution, or warm water</p> <p>Detach camera from scope (or camera from coupler), inspect and clean lens as needed</p> <p>Use suction or compressed air to dry out moisture (don't use cotton tip applicators on multi-pronged plug)</p> <p>Replace cables as necessary</p> <p>Move electrosurgical unit to different circuit or away from video equipment</p> <p>Reattach video cable at each monitor</p> <p>Adjust camera focus ring</p> <p>Inspect scope/camera, replace if needed</p> <p>Adjust enhancement and/or grain settings for units with this option</p>
<p>7. Inadequate suction/irrigation</p>	<p>Occlusion of tubing (kinking, blood clot, etc.)</p> <p>Occlusion of valves in suction/irrigator device</p> <p>Not attached to wall suction</p> <p>Irrigation fluid container not pressurized</p>	<p>Inspect full length of tubing. If necessary, detach from instrument and flush tubing with sterile saline</p> <p>Detach tubing, flush device with sterile saline</p> <p>Inspect and secure suction & wall source connector</p> <p>Inspect compressed gas source, connector, pressure dial setting</p>
<p>8. Absent or "weak" cauterization</p>	<p>Patient not grounded properly</p> <p>Connection between electro-surgical unit and instrument loose</p> <p>Foot pedal or hand switch not connected to electrosurgical unit</p> <p>Wrong output selected</p> <p>Connected to the wrong socket on the electrosurgical unit</p> <p>Instrument insulation failure outside of surgeon's view</p>	<p>Assure adequate grounding pad contact</p> <p>Inspect both connecting points</p> <p>Make connection</p> <p>Correct output choice</p> <p>Check that cable is attached to endoscopic socket</p> <p>Use new instrument and inspect insulation</p>