

# Face Protection Under N95 Masks

\*Used with permission from Stony Brook University

Product	Usage <b>T=</b> Treatment <b>P=</b> Prevention	Lawson#	Comfort	Ease of removal	Total Fit Factor*
Mepilex Light	<b>T</b> <b>P</b>	62444	Very	Easy no redness	148
Duoderm Extra Thin	<b>T</b> <b>P</b>	16684	Very	Easy Tenderness redness	167
Tegaderm	<b>P</b>	06527	Very Makes skin feel tight	Challenging Tenderness	186
Control: Cavilon *	<b>P</b>	01069	Very	Wears off over time	200

**\*Total Fit Factor: 100 or Greater = N95 is effective**

**\* Cavilon Skin Protectant Barrier recommended for use under each dressing type ( allow to completely dry before placing dressings)**



**Dressings should be placed using 3 pieces**

**Cheekbones: 1 inch x 2 inches**

**Bridge of Nose: 1 inch x 1 ½ inches**

**Dressings should sit under edges of the N95 Mask**